## 1 View Installation History

- The option in the LOTS Help menu called 'View Installation History' (Figure 1.1) displays a history of every installation performed on the Master and Slave PCs.
- 2. The most recent installations are shown at the top of the list (see Figure 1.2).

► Note: If you only want to show the installation history of the current PC, go to the 'View Installation History' on that PC and tick the '<u>This Computer</u> Only' tick box (as shown in Figure 1.2).

You can also filter for just 'Program Updates' or just 'Drug Updates' by using the tick boxes located at the top of the window (as shown in Figure 1.3).

 You can view the original PDF Release Notes by highlighting the installation in the list and clicking <sup>(</sup><u>V</u>iew Release Notes' (see Figure 1.3).

> Note: If no Release Notes are found for a selected release, the '<u>V</u>iew Release Notes' button will be disabled.

4. To print out this list click the '<u>P</u>rint' button.

🍯 LOTS Start Menu							
File	Help	LOTS					
	LO	TS Help	F1				
	Ver	rsion	Ctrl+V				
	Ab	out	Ctrl+A				
	View Installation History						
	Vie	w Outda	ited Modules				
Figur	e 1.1						

CLOTS Installation Histor	y			
🔽 Program Updates 🔽 🛛	<u>D</u> rug Updates	Double-click a column heading to sort the list. Double-click again to sort in reverse order.		
Installation Filename	Update Type	Installation Date	Computer Name	Master?
L2016R00.EXE	Program Update	01/04/2008 8:20 AM	TILL1	No
L2016R00.EXE	Program Update	01/04/2008 8:00 AM	DISPENSE1	Yes
L2015R00.EXE	Program Update	27/02/2008 8:03 AM	DISPENSE1	Yes
L2015R00.EXE	Program Update	27/02/2008 8:09 AM	TILL1	No
D0803R00.EXE	Drug Update	01/03/2008 8:12 AM	TILL1	No
D0803R00.EXE	Drug Update	01/03/2008 7:50 PM	DISPENSE1	Yes
☐ <u>T</u> his Computer Only		⊻iew Release Not	tes <u>P</u> rint	Close

Figure 1.2

✓ Program Updates I	<u>D</u> rug Updates	Double-click a column heading to sort the l Double-click again to sort in reverse or			
Installation Filename	Update Type	Installation Date	Computer Name	Maste	
L2019R00.EXE	Program Update	12/03/2009 10:27 AM	VICWKS006	Yes	
D0811R00.EXE	Drug Update	19/02/2009 11:06 AM	VICWKS006	Yes	
L2018R02.EXE	Program Update	23/01/2009 12:13 PM	VICLAP009	Yes	
L2018R00.EXE	Program Update	07/01/2009 11:16 AM	VICLAP009	Yes	
L2017R73.EXE	Program Update	10/12/2008 4:07 PM	VICLAP009	Yes	
L2017R70.EXE	Program Update	01/12/2008 9:47 AM	VICLAP009	Yes	
D0811R00.EXE	Drug Update	12/11/2008 12:32 PM	VICLAP009	Yes	
L2017R67.EXE	Program Update	11/11/2008 3:17 PM	VICLAP009	Yes	
L2017R18.EXE	Program Update	07/10/2008 9:46 AM	VICLAP009	Yes	
L2017R05.EXE	Program Update	19/09/2008 3:14 PM	VICLAP009	Yes	
L2017R00.EXE	Program Update	19/09/2008 11:33 AM	VICLAP009	Yes	
D0807R00.EXE	Drug Update	11/07/2008 10:01 AM	VICLAP009	Yes	
L2016R70.EXE	Program Update	11/07/2008 9:09 AM	VICLAP009	Yes	
L2016R65.EXE	Program Update	12/06/2008 4:09 PM	VICLAP009	Yes	
D0806R00.EXE	Drug Update	12/06/2008 3:52 PM	VICLAP009	Yes	
< 💷				>	

Figure 1.3